

ALVA *USB* 640

User Manual

Version 1.3

UDC-02067B

© 2010 Optelec, the Netherlands.
All rights reserved



Optelec
2993 LT Barendrecht
The Netherlands
Telephone: +31 (0)88 6783 444. Fax: +31 (0)88 6783 400
Internet: www.optelec.com

Contents

- 1. Introduction..... 2
 - 1.1. About this manual 2
 - 1.2. ALVA USB 640 package contents 2
- 2. Getting to know your ALVA USB 640 3
 - 2.1. ALVA USB 640 3
 - 2.1.1. Top and front panel layout 3
 - 2.1.2. Right side layout 4
 - 2.1.3. Left side layout..... 4
 - 2.1.4. Front panel keys, eTouch navigation buttons and cursor routing keys 4
 - eTouch keys 4
 - Cursor routing keys 4
 - Thumb keys 4
- 3. Installation 5
 - 3.1. Setting up the ALVA USB 640 5
 - 3.2. Installing the ALVA USB 640 screen reader drivers 5
- 4. ALVA USB 640 key assignments for JAWS, Dolphin and Window-Eyes 6
- Regulatory notice..... 8
- CE Caution 8
- Explanation of Graphical Symbols..... 9
- Safety Logo: 9
- Disposal of old products 9
- Product specifications..... 9
- FCC 10
- FCC RF Radiation Exposure Statement..... 10
- Optelec offices worldwide** 11

1. Introduction

Congratulations on choosing the ALVA USB 640, a tactile user interface with high quality Braille cells. The ALVA USB 640 offers optimal flexibility and productivity in an ultra compact design. Used in combination with screen reading software, the ALVA USB 640 can be installed for specific operating systems.

The ALVA USB 640 extends the navigational and operational possibilities of your screen reading software, and allows you to be in full control of your computer. Its ergonomical design allows for optimal efficiency when working with a computer. Note that the ALVA USB 640's functionality will depend on the Braille display options of your screen reader and your mainstream computing hardware.

The ALVA USB 640 was designed to adhere to the Open Standard philosophy so that it may be used with mainstream technology, enabling compatibility with most screen-reading software and operating systems. For the latest information about screen readers that fully support the device, or for any questions or suggestions concerning the use of this Braille display, please contact your supplier or the Optelec office nearest you by consulting the contact information page at the end of this manual. Your feedback is highly appreciated. We hope you enjoy working with your ALVA USB 640.

1.1. About this manual

This manual will familiarize you with the features and operation of your ALVA USB 640. Please read this manual carefully before using your device.

1.2. ALVA USB 640 package contents

The ALVA USB 640 package contains the following:

- The ALVA USB 640 Base Unit
- Long and short USB cables
- The ALVA BC driver software for Windows XP, Windows Vista and Windows 7, and product documentation on CD-ROM.

If any of these items are missing from your package, please contact your supplier.


2. Getting to know your ALVA USB 640

This chapter will provide you with a detailed description of the ALVA USB 640 layout.

2.1. ALVA USB 640

The ALVA USB 640 features:

1. 40-cell Braille display
2. eTouch keys (left) for panning left and right, and for screen-reading functions
3. eTouch keys (right) for panning left and right, and for screen-reading functions
4. Thumb keys for panning left, right, up and down, and for returning to the cursor position
5. 40 Cursor Routing Keys
6. Mini USB connector
7. On/off switch



When the ALVA USB 640 is placed correctly in front of you, the Braille cells are located at the back of the flat surface, and the thumb keys are facing you.

2.1.1. Top and front panel layout

The ALVA USB 640 design features a row of 40 Braille cells, located at the back of the flat surface. In front of the Braille cells, you will find a row of small buttons, which are known as cursor routing keys. Each Braille cell has a corresponding cursor routing key that can be used for multiple functions. On the left and right side of the Braille line, you will find two groups of two keys. These are called eTouch keys, which can be used for Braille panning left and right, and also for various other screen-reading functions.

2.1.2. Right side layout

The power switch is located on the right side of the ALVA USB 640. The ALVA USB 640 is switched on when the rectangular push-button is pressed. Upon switching it on, the ALVA USB 640 will generate an audio signal.

2.1.3. Left side layout

The left side of the ALVA USB 640 features two connectors. From front to back, you will find a mini USB connector and a feature connector. The feature connector has no function with the release of this manual.

2.1.4. Front panel keys, eTouch navigation buttons and cursor routing keys

The names of the ALVA USB 640 keys listed in this section will be used throughout all of the documentation, including the screen reader documentation. This section will present the name of the keys, and describe their main function as supported by the screen reader. Their precise function is defined by the screen reader. Please refer to your screen reader documentation for detailed key function definitions.

eTouch keys

The two keys found on each side of the Braille line are called eTouch keys. They are used for panning left and right, and also for various other Braille functions.

- eTouch 1 Top left key marked with two vertical lines
- eTouch 2 Bottom left key marked with one horizontal line
- eTouch 3 Top right key marked with two vertical lines
- eTouch 4 Bottom right key marked with one horizontal line

Cursor routing keys

Each cursor routing key corresponds to the Braille cell located immediately above it. These keys are typically used for routing the cursor or caret to the corresponding Braille cell, generating a left or right mouse click at that point, or for obtaining additional information about the content of the Braille cell.

Thumb keys

At the front of the ALVA USB 640, you will find five front panel keys known as thumb keys. These include four large keys symmetrically aligned along the front of the ALVA USB 640, and one smaller button located in the centre. Pressing the four large thumb keys will allow you to navigate left, right, up or down respectively. Pressing the middle key will allow you to return to the cursor position.

3. Installation

This section will present a general description of the installation procedure for the JAWS, Dolphin and Window-Eyes screen readers on a Windows XP, Windows Vista or Windows 7 platform. To install and configure the Braille display with other screen readers and operating systems, please refer to the screen reader documentation for instructions.

3.1. *Setting up the ALVA USB 640*

The ALVA USB 640 is compatible with Windows XP, Windows Vista and Windows 7. Previous versions of Windows are not supported. Windows XP, Windows Vista or Windows 7 will automatically detect the ALVA USB 640, and install the necessary USB drivers.

1. Attach the USB cable to the ALVA USB 640 and then to the computer.
2. Switch on the ALVA USB 640. Windows will automatically install the necessary USB drivers. This does not require any further user input. Once the connection is established, the ALVA USB 640 will generate an audible signal, and will display the message "USB connected".

3.2. *Installing the ALVA USB 640 screen reader drivers*

1. Insert the ALVA USB 640 software CD-ROM that came with the package in your computer's disk drive.
2. The program will automatically run if the autostart feature of your CD-ROM drive is enabled. If the autostart feature is not enabled, open the Windows Start Menu; select "Run" and type d:\setup.exe (assuming d: is the letter assigned to your disk drive); and then press Enter.
3. Select "Install screen reader drivers".
4. The setup program will guide you through the installation.

You have the choice between two setup options: Complete or Custom setup. We recommend that you use the Complete setup option.

5. a) The Complete setup will automatically install and update all the necessary components.
5. b) If you would like to install the ALVA USB 640 software manually, you may choose the Custom setup option. This option is for advanced users only. The Custom setup will ask you what components you would like to install. For each component chosen, you will need to provide the correct installation location.

Drivers for screen readers other than JAWS, Dolphin and Window-Eyes are included in the products themselves, and are not found on the ALVA USB 640 CD-ROM. Please refer to your screen reader documentation for more information about the ALVA USB 640 setup.

4. ALVA USB 640 key assignments for JAWS, Dolphin and Window-Eyes

T = Thumb key

ALVA USB 640 key	JAWS	Dolphin Supernova / Hal
eTouch1	Left wiz wheel up	Braille document read back
eTouch2	Left wiz wheel down	Where am I
eTouch3	Right wiz wheel up	Braille document read forward
eTouch4	Right wiz wheel down	Read status line
eTouch1 + eTouch2	Press left wiz wheel	Braille line start
eTouch1 + eTouch3	Braille auto advance	
eTouch1 + eTouch4	Increase auto advance speed	
eTouch2 + eTouch3	Decrease auto advance speed	
eTouch2 + eTouch4	Braille toggle speech interrupt	
eTouch3 + eTouch 4	Press right wiz wheel	Braille line end
T3 + eTouch1	Braille toggle mode	Physical mode on/off
T3 + eTouch2	Braille toggle Cursor	Toggle virtual focus on / off
T3 + eTouch3	Braille moves active	Braille verbosity scheme cycle
T3 + eTouch4	Braille toggle 8 pixels per space	Expand gaps on/off
T1	Pan left	Braille back
T2	Braille prior line	Line up
T3	Route Braille to active cursor	Go to focus
T4	Braille next line	Line down
T5	Pan right	Braille forward
T1 + T2	Braille Top	Braille top
T1 + T3	Set Braille verbosity	
T1 + T4	Braille change cursor shape	Cursor style
T1 + T5	Braille toggle table reading	Braille input on/off
T2 + T4	Braille toggle 8 dots	Toggle 8 dots/6 dots
T2 + T3	Braille Grade2 translation	Literary Braille on/off
T2 + T5	Braille toggle marking off	Braille cursor visible on/off
T3 + T4	Braille Grade2 expand current word	
T3 + T5	Braille toggle marking	Enhance on/off
T4 + T5	Braille bottom	Braille bottom

ALVA USB 640 key	Window-Eyes
eTouch1	Scroll Braille left
eTouch2	Top line
eTouch3	Scroll Braille right
eTouch4	Bottom line
T3 + eTouch1	Scroll mode toggle
T3 + eTouch2	WE/Mouse toggle
T3 + eTouch3	Auto route cursor toggle
T3 + eTouch4	Spacing toggle
T1	Scroll Braille left
T2	Prior Braille line
T3	To focus
T4	Next Braille line
T5	Scroll Braille right
T1 + T2	Show attribute toggle
T1 + T3	Auto route cursor toggle
T1 + T4	Indent toggle
T1 + T5	Quick message rotor
T2 + T4	6 / 8 dot toggle
T2 + T3	Grade 2 Braille toggle
T2 + T5	Braille tracking toggle
T3 + T4	Auto untranslate at cursor
T4 + T5	Attribute to show rotor

Regulatory notice

This equipment has been developed and tested in order to comply with the following regulations:



and  for a Class B digital device, pursuant to part 15 rules.

Medical electrical equipment requires special precaution regarding EMC (Electro Magnetic Compatibility), and must be installed and put into service according to EMC guidelines.

WARNING: The use of accessories, transducers and cables other than those specified and sold by the manufacturer of the ALVA USB 640 as replacement component parts may result in increased emissions or decreased immunity of the ALVA USB 640.

By design, the ALVA will shut down in case of an ESD 8kV test pulse. The ALVA USB 640 will have to be switched back on to resume operation.

CE Caution

- The input and output parts need to be connected correctly;
- The unit is designed for exclusive interconnection with IEC/EN 60XXX certified equipment outside of the user environment, and with IEC/EN 60601-1 certified equipment within the user environment.
- This device complies with EN 60601-1-2. To minimize interference from other equipment, a minimum distance of 0.5 m shall be kept from other potential electromagnetic sources, such as a mobile phone.
- Equipment connected to analog or digital interfaces must comply with the respective IEC/EN standards (e.g. IEC/EN 60950-1 for data-processing equipment, and with IEC/EN 60601-1 for medical equipment).
- Do not open the housing.
- Only qualified service personnel should service the unit.
- All possible forms of configuration must comply with the current version of the IEC/EN 60601-1-1 system standard directive.
- Please be aware that connecting additional equipment to the signal input part or signal output part of the device entails that a medical system is being configured. Therefore, the individual, who is configuring and interconnecting the system, is responsible for the entire system, and must make sure it is compliant with the current version of the requirements of the IEC/EN 60601-1-1 system standard. In case of any doubt about system compliance under the current version of the requirements of the IEC/EN 60601-1-1 system standard, please consult the technical service department or your local representative.

Explanation of Graphical Symbols



This symbol warns the user that uninsulated voltage within the unit may have sufficient intensity to cause electric shock, therefore it is dangerous to make contact with any part inside this unit.



This symbol alerts the user that important literature concerning the operation and maintenance of this unit is included. It should therefore be read carefully in order to avoid any problems.



Stand-by



Direct Current

Safety Logo:



This device complies with the 93/42/EEC, EN60601-1, EN60601-2 relative to European standards.

Disposal of old products

-If you see this symbol-

Information on Disposal methods for other countries outside the European Union.



This symbol is only valid in the European Union.

If you wish to discard this product, please contact your local authority or dealer to learn about the correct disposal method.

Product specifications

Specifications

Environment	Operating temperature: +0° C ~ +40° C
Conditions	Storage/transportation temperature: -20° C ~ +60° C
	Operating relative humidity: 20% ~80% non-condensing
	Storage/transportation relative humidity: 10% ~95% non-condensing
	Atmospheric pressure: 500~1060 hPa

FCC

FCC Notice

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and the receiver.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Caution

To ensure continued compliance, follow the attached installation instructions, and use only shielded interface cables with ferrite core when connecting to a computer or peripheral device. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

FCC RF Radiation Exposure Statement

This transmitter must not be co-located or operated in conjunction with any other antenna or transmitter. This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. To maintain compliance with FCC RF exposure compliance requirements, avoid direct contact to the transmitting antenna during transmitting. Any changes or modifications made to this device (including the antennas) not expressly approved by the manufacturer may void the user's authority to operate the equipment.

Warranty conditions

Optelec warrants the ALVA USB 640, effective from the date of delivery, to be free of any defects in material and workmanship.

Warranty is not transferable and does not apply to groups, multi users or agencies.

The ALVA USB 640 has been designed for the individual purchaser to be used in home and office environments. Optelec reserves the right to repair or replace any purchased ALVA USB 640 with a similar or better product.

In no event shall Optelec or its distributors be liable for indirect or consequential damages. The original user's remedies are limited to replacement of ALVA USB 640 modules. This warranty is valid only when serviced in the country of original purchase, and with intact seals. For additional warranty claims or service during or after the warranty period, please consult your supplier.

Optelec takes no responsibility for use of this device other than described in this manual.

Optelec offices worldwide

Optelec U.S. Ltd

3030 Enterprise Court, Suite C
Vista, CA 92081
U.S.A.
T. +1 (800) 826-4200
www.optelec.com

Optelec Canada

1832, Boulevard Marie-Victorin
Longueuil, QC J4G 1Y9
Canada
T: +1 (800) 665-3005
www.optelec.com

Optelec Limited

2, Millfield House, Woodshots Meadow
Croxley Green Business Park
Watford, WD18 8YX
United Kingdom
T: +44 (0)1923 23 13 13
www.optelec.co.uk

Optelec Nederland BV

Breslau 4
2993 LT Barendrecht
The Netherlands
T: +31 (0)88 6783 555
www.optelec.nl

Optelec NV

Baron Ruzettelaan 29
B-8310 Bruges
Belgium
T: +32 (0)50 35 7555
www.optelec.be

Optelec GmbH

Fritzlarer Strasse 25
34613 Schwalmstadt
Germany
T: +49 (0)6691 96170
www.optelec.de

Optelec Corporate Headquarters: Breslau 4 - 2993 LT Barendrecht – The Netherlands – T: +31 (0)88 678 3444 – www.optelec.com

www.optelec.com

