

Copyright © 2007 Optelec B.V.

Breslau 4

2993 LT Barendrecht

The Netherlands

Phone: +31 88 678 34 44

Fax: +31 88 678 34 00

E-mail: [dealersupport@optelec.nl](mailto:dealersupport@optelec.nl)

Website: <http://www.optelec.com>

## Preface

Congratulations on choosing the ALVA 584 Satellite Pro Braille display. The 584 Satellite Pro is the fifth-generation Braille display, featuring the latest hardware innovations and ergonomic designs. Equipped with inbound and outbound USB ports, an on-board multi-voltage power supply, Touch Cursor and Double Touch Cursor buttons, and a redesigned front panel, the ergonomic design of the ALVA 584 Satellite Pro gives the user extended flexibility and productivity, all housed in a new, thin, light-weight profile.

This manual will familiarize you with the features and basic operation of the ALVA 584 Satellite Pro. Please read this manual carefully before connecting the Braille display to your computer.

If you have any questions or suggestions concerning the use of this Braille display, please contact your distributor or Optelec's home office using the contact information below. Your feedback is greatly appreciated. We hope you enjoy working with your ALVA 584 Satellite Pro.

Optelec B.V.

Breslau 4

2993 LT Barendrecht

The Netherlands

Phone: +31 88 678 34 44

Fax: +31 88 678 34 00

E-mail: [dealersupport@optelec.nl](mailto:dealersupport@optelec.nl)

Website: <http://www.optelec.com>

## Regulatory notice

This equipment has been developed and tested in order to comply with the following regulations:

- CE
- FCC for a Class B digital device, pursuant to part 15 rules

These limits are designed to provide reasonable protection against harmful interference in a residual installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Re-orient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the distributor or an experienced radio/TV technician for help.

# Contents

<b>Preface</b>	<b>iii</b>
<b>Regulatory notice</b>	<b>iv</b>
<b>Contents</b>	<b>v</b>
<b>1 Introduction</b>	<b>7</b>
1.1 About this manual	7
1.2 ALVA 584 Satellite Pro package contents	8
<b>2 Description of the ALVA 584 Satellite Pro</b>	<b>8</b>
2.1 Surface and front panel layout	8
2.2 Back panel layout	8
2.3 Front panel keys, satellite keypads and double touch cursors	8
2.3.1 Front panel keys	8
2.3.2 Windows satellite keypad	8
2.3.3 Screen reader satellite keypad	8
2.3.4 Touch cursors and double touch cursors	8
<b>3 Installation</b>	<b>8</b>
3.1 USB connection	8
3.1.1 Installing the USB driver for Microsoft Windows	8
3.2 Serial connection	8
3.2.1 Serial port	8
3.2.2 Serial cable	8
3.2.3 Installing the serial driver for Microsoft Windows	8
3.3 Switching between USB and serial connection	8
<b>4 Local mode menu</b>	<b>8</b>
4.1 Activating the Local mode menu	8
4.2 Main menu	8

4.3	Status submenu	8
4.4	Config submenu	8
4.5	Save submenu	8
4.6	Load submenu	8
4.7	Test submenu	8
<b>5</b>	<b>Miscellaneous</b>	<b>8</b>
5.1	Error Messages	8
5.2	Using both the serial and USB connections	8
<b>Index</b>	<b>8</b>	

# 1 Introduction

The ALVA Satellite Braille display presents text, graphic and other information in Braille corresponding to what is displayed on a computer screen. Used together with screen reader software the Satellite display can be installed for a particular operating system such as Windows 98, Windows 2000 or Windows XP to provide a tactile and audible user interface.

The ALVA 584 Satellite Pro Braille display offers an enhanced ergonomic design that extends the navigation and operation possibilities of the screen reader software. The features of the ALVA 584 Satellite Pro are designed for optimal capabilities when working in a graphical computer environment. Note that the functionality of the ALVA 584 Satellite Pro depends somewhat on the Braille display options that your screen reader software offers.

The Satellite Braille display was designed to adhere to Open Standard philosophy. This enables compatibility with most screen reader software and operating systems. Contact your distributor or Optelec for updated information about screen reader software that fully supports the functionality and possibilities of the ALVA 584 Satellite Pro Braille display or tailor-made solutions.

ALVA Satellite Braille displays are optimized to work in a graphical computer environment. It is also possible to work in a non-graphical environment under MS-DOS.

## 1.1 About this manual

This manual is intentionally kept brief to enable you to quickly begin using your ALVA Satellite Braille display and for easy reference in the future. Chapter 1 provides an introduction and a list of the contents

of the ALVA 584 Satellite Pro package. Chapter 2 provides a detailed description of the ALVA 584 Satellite Pro and the names used for the keys, buttons and switches on the Braille display. Chapter 3 describes how to set up and install the Braille display. Chapter 4 explains how to configure the Braille display. Chapter 5 provides information regarding the integrated power supply. Chapter 6 provides information regarding the Universal Serial Bus connections and explains how to operate the ALVA 584 Satellite Pro in an MS-DOS environment. If you need additional information about the ALVA 584 Satellite Pro Braille display, please consult your distributor.

## 1.2 ALVA 584 Satellite Pro package contents

The ALVA 584 Satellite Pro package contains the following:

**Box 1:** – The ALVA 584 Satellite Pro Braille display

**Box 2:** – An AC-adapter  
– Serial cable to connect the Braille display to a computer  
– USB cable to connect the Braille display to a computer  
– The ALVA Satellite driver software on CD-ROM  
– This manual in print and Braille

If any of these items are missing from your package, please contact your distributor.

## 2 Description of the ALVA 584 Satellite Pro

### 2.1 Surface and front panel layout

The ALVA 584 Satellite Pro Braille display is designed with a row of 84 Braille cells at the front of the flat surface. Behind the Braille cells, a flat surface is provided so that a keyboard can be placed on top of the Braille display. Above the row of Braille cells are two rows of small buttons known as Touch Cursors (lower row) and Double Touch Cursors (upper row). Each Braille cell has a corresponding Touch Cursor and Double Touch Cursor button.

On the left and right side, just below the array of Braille cells, are two groups of six keys. These groups are called the Satellite keypads. Each Satellite keypad forms an inner circle of four keys with two larger keys positioned alongside the circle. At the front of the Braille display are ten front panel keys: two round keys, one on the left and one on the right, four rectangular keys in the middle, and aside two tumbler keys.

### 2.2 Back panel layout

At the back of the Braille display are the power switch and connections for the PC, other devices and the power supply. The connections are partially covered by the case, which extends somewhat beyond the back panel. Tactile markings have been placed on the underside of the cover that indicate the meaning of the connections.

From left to right (with the Braille display front panel facing you) the layout of the connections is:

- On/off switch, indicated by a circle and a dash. The Braille display is off when tumbled to the left.
- One 9-pins serial connector for serial connection to the PC, indicated by two circles between vertical marks
- Inbound USB connector, square shaped, for connection to the PC, indicated by an inward pointing arrow to the connector
- Electrical socket connector indicated by a curl

## 2.3 Front panel keys, satellite keypads and double touch cursors

This section lists the standard names for all front panel keys and satellite keypads. These names will also be used in the user manual of your screen reading software when referring to the ALVA 584 Satellite Pro. The exact functionality of the keys is defined by the screen reader. Please refer to your screen reader documentation for detailed functionality of the ALVA 584 Satellite Pro.

### 2.3.1 Front panel keys

From left to right, the front panel keys are:

Left tumbler key

HOME Round front panel key at the left

LEFT Rectangular key with bend at the left

UP Rectangular key with line mark at the top

DOWN Rectangular key with line mark at the bottom

RIGHT Rectangular key with bend at the right

CURSOR Round front panel key at the right

Right tumbler key

### **2.3.2 Windows satellite keypad**

One of the two satellite keypads will be used for issuing Windows commands. By default this will be the satellite keypad on the right hand side of the Braille display. The following names will be used in the documentation:

BACK TAB	Left outer circular key
FORWARD TAB	Right outer circular key
LEFT PAD	Left inner circular key
UP PAD	Top inner circular key
DOWN PAD	Bottom inner circular key
RIGHT PAD	Right inner circular key

### **2.3.3 Screen reader satellite keypad**

The second keypad is defined by the functionality of the screen reader. Examples of screen reader functionality that can be implemented on this satellite keypad are speech-related functions, extended Braille functionality, and toggles for user defined scripts and macros. Your screen reader documentation describes the functionality that is implemented in the screen reader Satellite keypad. By default the satellite keypad located at the left-hand side of the ALVA 584 Satellite Pro will be assigned to the screen reader functionality.

The following names will be used in the documentation:

BACK EAR	Left outer circular key
FORWARD EAR	Right outer circular key
WEST PAD	Left inner circular key
NORTH PAD	Top inner circular key

SOUTH PAD      Bottom inner circular key

EAST PAD        Right inner circular key

#### **2.3.4 Touch cursors and double touch cursors**

The touch cursors and double touch cursors correspond to the particular Braille cell immediately below each pair of buttons. Typically, they are used for routing the cursor or caret to the corresponding Braille cell, generating a left or right mouse click at that point, or for obtaining additional information about the contents of the Braille cell.

The touch cursors are referred to as:

TC 1-84 This is the row of 84 buttons directly above the Braille cells

The double touch cursors are referred to as:

DTC 1-84 This is the top row of 84 buttons above the Braille cells

## 3 Installation

In this section, the installation procedure is only described in general. The file 'ReadMe' (located at the root of the ALVA Braille Display Drivers CD-ROM included in the ALVA Satellite Pro package) contains a detailed description of the installation procedure. The 'ReadMe' file describes any deviations from the procedure explained below and also describes the installation for some of the popular screen readers.

### 3.1 USB connection

#### 3.1.1 Installing the USB driver for Microsoft Windows

Before starting the installation process, make sure the ALVA Satellite is switched off. Upon Windows start-up, insert the ALVA Braille Display Drivers CD-ROM that came shipped with the ALVA Satellite package. When you have connected the Satellite to the PC using the USB cable, Windows will detect the integrated Generic USB Hub (see also section **Error! Reference source not found.**), start the 'Add New Hardware Wizard' and automatically install the appropriate driver. This does not require any user input. After having inserted the CD-ROM, the installation program will automatically run if the autostart feature for your CD-Rom drive is enabled. If not, open 'Run' from the start menu, type d:\setup (assuming d: is the letter assigned to your CD-ROM drive) and press Enter. The installation program will prompt you for the folder to where it should copy the drivers. Typically this is the main folder of your screen reader (e.g. in the case of JAWS for Windows) or one of its subfolders (e.g. in the case of Supernova). Specify the path and press Enter. When the installation program has finished, it will ask you to restart your

system. Press Enter to do so, or select "No, I will restart my computer later" and press Enter to return to Windows. Note that the installation will not be completed until the system has been restarted.

When switching on the Satellite after having restarted the PC, the Windows Add New Hardware Wizard will appear and automatically complete the installation. This does not require any further user input. The ALVA Satellite will then issue an audible signal to indicate that the connection is working.

## 3.2 Serial connection

### 3.2.1 Serial port

Use the serial port for communication with operating systems that do not support USB (MS-DOS, Windows 95 or Windows NT4). The Braille display automatically switches to serial communication when the USB connection is not used.

### 3.2.2 Serial cable

Use one of the serial cables supplied with the ALVA 584 Satellite Pro Braille display for connection to the computer. Plug the serial cable in to the 9-pin serial port of the Braille display (indicated by two circles between vertical marks) and connect the other end of the cable to a serial port of the PC.

### 3.2.3 Installing the serial driver for Microsoft Windows

Upon Windows start-up, insert the CD-ROM that came shipped with the ALVA Satellite package. The installation program will automatically run if the autostart feature for your CD-Rom drive is enabled. If not, open 'Run' from the start menu, type d:\setup (assuming d: is the letter assigned to your CD-ROM drive) and press Enter. The installation program will then prompt you for the folder to

where it should copy the drivers. Typically this is the main folder of your screen reader (e.g. in the case of JAWS for Windows) or one of its subfolders (e.g. in the case of Supernova). Specify the path and press Enter. When the installation program has finished, it will ask you to restart your system. Press Enter to do so, or select "No, I will restart my computer later" and press Enter to return to Windows. Note that the installation will not be completed until the system has been restarted.

### 3.3 Switching between USB and serial connection

The ALVA 584 Satellite Pro automatically detects whether it is connected to the serial port or the USB port of a computer. The ALVA 584 Satellite Pro can be connected to two different computers simultaneously. One computer will use the USB port and the other the serial port. In this case the ALVA 584 Satellite Pro will choose the USB port as the main communication port. In the Local mode menu of the ALVA 584 Satellite Pro you can specify whether the USB or serial port is preferred.

Open the Local mode menu by simultaneously pressing the Back Ear, South Pad and Forward Ear keys of the Screen reader Satellite keypad. The Braille display shows the text *usb ser* and the other menu options. Press the Touch Cursor above the *usb* or *ser* option to select the preferred connection. This connection will be used next time you power on the ALVA 584 Satellite Pro.

Exit the Local mode menu by pressing the CURSOR front panel key.

## 4 Local mode menu

When the ALVA Satellite is turned on, the following message is displayed:

ALVA Satellite OK, waiting for connection... This message disappears when the screen reader is loaded and a successful connection to the Satellite is established.

### 4.1 Activating the Local mode menu

The ALVA Satellite menu is called the Local mode menu, and it enables you to change several built-in settings of the Braille display. The Local mode menu can be activated by pressing either of the following key combinations: Back Ear + South Pad + Forward Ear of the Screen reader Satellite keypad or Back Tab + Down Pad + Forward Tab of the Windows Satellite keypad.

When the Local mode menu is activated, a sound is emitted and the main menu appears on the Braille display:

```
usb ser status config save load test
```

The relevant submenu is activated by pressing one of the Touch Cursors above the main menu option. Use the HOME front panel key to move back to the previous menu level or to exit the Local mode menu when the main menu is active.

Use the CURSOR front panel key to immediately exit the Local mode menu from anywhere within the Local mode.

## 4.2 Main menu

```
usb ser status config save load test
```

The main menu has seven options:

- usb For manually setting the USB connection. Select this option when the USB port must be used for communication with the PC.
- ser For manually setting the serial connection. Select this option when the serial port must be used for communication with the PC.
- status Contains a submenu with the current settings: connection, battery information, serial number and firmware date of the ALVA 584 Satellite Pro (see section 4.3).
- config Contains a submenu for changing a variety of settings on the ALVA Satellite (see section 4.4).
- save Contains a submenu for saving the current settings (see section 4.5).
- load Contains a submenu for loading particular settings or restoring the default settings (see section 4.6).
- test Activates a test mode used for testing the operation of the ALVA Satellite Braille display (see section 4.7).

## 4.3 Status submenu

The status submenu is used to verify the current status of the ALVA Satellite. The status menu contains the menu items:

```
comm battery sernum settings
```

**comm** Refers to the status of the communication. Activating the comm menu item provides the status of the USB or serial communication, and shows the device that is attached to the serial connection. This can be an *ALVA spacepad*, *keypad* or *FeelMe utility*, but is generally a *PC*. The most frequently displayed communication information is:

```
usb mode=pc
```

```
ser mode=pc
```

**battery** Does not apply for the ALVA 584 Satellite Pro. The battery status shows:

```
batteries not inserted
```

**sernum** The serial number of the Braille display and the version number of the firmware are displayed. The firmware is loaded into the main processor of the ALVA 584 Satellite Pro and is responsible for the functionality of the Braille display. If technical developments require new functionality, the current firmware version may need to be replaced by a newer version. An example of the sernum status is:

```
serial nr=xxxxxxxx firmware version=2.1.0.0
```

where xxxxxxxx stands for the serial number, such as 58400001.

**settings** The settings menu item provides the status of the current configuration of the ALVA satellite. Use the UP and DOWN front panel keys to view the settings. The settings menu is used only for viewing the current settings. If you want to change one of these configuration settings use the config menu (section 4.4). The settings menu provides the following information:

```
loaded set=xxxx
```

where xxxx indicates which set, set1 to set7, is currently loaded.

```
volume=xx
```

where xx is a number from 0 (no sound) to 31 (maximum volume) referring to the volume of the built-in speaker.

```
sleep mode=xxx
```

The sleep mode status can be set from 1 second to 255 seconds. This indicates the time after which the Braille display puts itself into sleep mode if there is no activity. Consequently, when the Satellite is in sleep mode, it consumes virtually no energy. The sleep mode is inactive if 0 is shown.

```
dot pressure=x
```

The dot pressure can be set from 1 to 5 indicating the amount of pressure that the dots produce when reading the Braille. A value of 1 indicates the least amount of dot pressure and 5 indicates the highest dot pressure.

```
status cells=xxxxx
```

The position of the three status cells can be left, right or none. This indicates whether three status cells and one blank separator cell are positioned on the left, right, or turned off. If *none* is selected, then all 84 Braille cells are used as data cells.

```
switch satellite function=xxx
```

This setting can be *yes* or *no* indicating whether the Windows and Screen reader satellite keypads are interchanged or not.

```
front panel=xxxxxxxx
```

This option displays the layout of the UP, DOWN, LEFT and RIGHT front panel keys, and also indicates whether or not the PROG key is emulated on one of the satellite keypads for compatibility with the ALVA DOS software.

```
internal Braille table=xxxxxx
```

This indicates the internal Braille table language that is used in the Local mode menu and the built-in messages that can appear on the Braille display.

## 4.4 Config submenu

The configuration submenu has seven user-definable options:

```
vol sleep dots stcell satellite front table
```

Select the option that you want to change by pressing the Touch Cursor above the menu item. Use the UP and DOWN front panel keys to change the setting. Changes to the settings for the loaded configuration set are immediately saved when you press the HOME or CURSOR front panel key to go back to the previous menu level or to exit the Local mode menu, respectively.

```
volume=xx
```

**vol** Changes the volume level of the sounds generated by the ALVA 584 Satellite Pro in steps from 0 to 31. The DOWN front panel key increases the volume; the UP front panel key decreases the volume. If you want no sounds to be generated, select 0. The default value is 16.

```
sleep mode=xxx
```

**sleep** Set the time from 1 second to 255 seconds, indicating the time after which the Braille display puts itself into sleep

mode if there is no activity. Consequently, the Satellite consumes virtually no energy when the Braille display is in sleep mode. The sleep mode is inactive when 0 is selected. The DOWN front panel key increases the time; the UP front panel key decreases the time. The default value is 255.

```
dot pressure=x
```

**dots** The dot pressure can be set from 1 to 5 indicating the amount of pressure that the dots provide when reading the Braille. The DOWN front panel key increases the dot pressure; the UP front panel key decreases the dot pressure.

```
status cells=xxxxx
```

**stcells** Select the position of the status cells with the UP and DOWN front panel keys. The options are *Left*, *Right* or *None*. *None* means that there are no status cells, and all Braille cells are used for displaying data.

```
switch satellite functions=xxx
```

**satellite** Interchange the Windows and Screen reader satellite keypads. The default setting is *no*, meaning that the Windows satellite keypad is at the right-hand side and the Screen Reader satellite keypad is at the left-hand side. Select *yes* if you prefer that the functionality be interchanged. We recommend that you choose a setting that leaves your preferred hand for Braille reading on the

Braille cells as much as possible so that the Windows Satellite keypad can be operated with your other hand.

```
front panel=xxxxxxxx
```

front Select the layout of the front panel keys. There are four options:

*LEFT UP DOWN RIGHT*, which is the default

*UP LEFT RIGHT DOWN*, which is the classic ABT3 and Delphi layout

*LEFT UP DOWN RIGHT* and *PROG*

*UP LEFT RIGHT DOWN* and *PROG*

which are the settings used to emulate the PROG key on the BACK Ear, enabling the ALVA Satellite to be fully compatible with ALVA DOS software.

```
internal Braille table=xxxxx
```

table Select the Braille table you want to use for the Local mode menus. Available options are *German*, *USA*, *German lc* (lower case) and *USA lc* (lower case). The first two Braille tables are 8-dot Braille and the last two are 6-dot Braille tables. Note that due to the limited memory of the processor, it is not possible to change the language of the Local mode. The Local mode menu is only available in English.

## 4.5 Save submenu

The configuration can be saved to seven different sets. In situations where several people use the same ALVA 584 Satellite Pro, it may be convenient that your particular configuration is saved to be reloaded when needed. The save submenu contains the options *set1* to *set7*.

```
set1 set2 set3 set4 set5 set6 set7
```

Press the Touch Cursor above the set in which you want to save the configuration. The changes are always saved automatically for the currently loaded set when you exit the Local mode menu.

## 4.6 Load submenu

The settings can be loaded by pressing the Touch Cursor above one of the seven sets. The last option in the menu is *def*, in which the predefined default configuration is restored under the name of the currently loaded set. The other sets will not be restored as the default.

```
set1 set2 set3 set4 set5 set6 set7 def
```

Press the Touch Cursor above the set of configurations you want to load. The loaded set can be viewed in the status submenu - option settings.

## 4.7 Test submenu

Use the test submenu if you have the impression that certain keys or Braille cells are not responding properly. The ALVA 584 Satellite Pro will activate a test mode for checking the Braille cells and function of

the Braille display keys. Pressing a particular key provides a unique dot pattern on the array of Braille cells:

Touch Cursor	Dots 3678 are shown at the cell location
Double Touch Cursor	Dots 1245 are shown at the cell location
Back ear	Dot 1 and 4 are shown in all cells
North pad	Dot 1 is shown in all cells
East pad	Dot 2 is shown in all cells
South pad	Dot 3 is shown in all cells
West pad	Dot 7 is shown in all cells
Forward ear	Dot 2 and 5 are shown in all cells
Back Tab	Dot 3 and 6 are shown in all cells
Up pad	Dot 4 is shown in all cells
Left pad	Dot 5 is shown in all cells
Down pad	Dot 6 is shown in all cells
Right pad	Dot 8 is shown in all cells
Forward Tab	Dot 7 and 8 are shown in all cells
Left front panel key	Dots in rows 1, 2, 3 and 7 are shown
Up front panel key	Dots in rows 4, 5, 6 and 8 are shown
Down front panel key	All dots are shown
Right front panel key	Indicates the cells for which all Touch Cursors and Double Touch Cursors were pressed while in test mode.
Left tumbler key to left	All dots odd cells are shown
Left tumbler key to right	All dots even cells are shown
Home	Back to previous menu in local mode

Cursor	Leave local mode and go back to host mode
Right tumbler key to left	Dots 1245 of even cells and dots 3678 of odd cells are shown
Right tumbler key to right	Dots 1245 of odd cells and dots 3678 of even cells are shown

## 5 Miscellaneous

### 5.1 Error Messages

When switching on the ALVA Satellite the Braille display tests whether various components of the internal hardware are working correctly. When the ALVA Satellite tests successfully, the following message is shown on the Braille display:

```
ALVA Satellite OK, waiting for connection...
```

This message disappears when a screen reader is loaded.

Various error messages can be presented on the Braille display. Pressing any key on the Braille display clears the error message and the ALVA Satellite may start depending on the nature of the message. These messages should never occur when the ALVA Satellite is properly functioning. Always consult your distributor when one of the following messages appears more than once:

```
Net voltage out of range
```

This indicates that there is a problem with the built-in power supply.

```
USB controller supply voltage out of range
```

The voltage of the USB controller cannot be set correctly, which indicates that the USB controller is not working properly.

Temperature sensor defect

The temperature sensor is defective which protects the ALVA Satellite against heating.

Serial eeprom defect

The eeprom in which the parameters of the ALVA Satellite are stored is defect.

Key shift register defect

One of the front panel keys, satellite keys or Double Touch Cursor prints is not working.

## 5.2 Using both the serial and USB connections

The ALVA Satellite can be connected with either the serial or USB connection to the PC, or with the USB connection to one PC and the serial connection to a second PC. The ALVA Satellite will automatically use the connection that is used for communication at that moment. For example, if the Braille display detects serial communication in a DOS environment, the serial connection will be automatically selected.

However, if both serial and USB communication is set up (for example when both PCs send data to the Braille display) the ALVA 584 Satellite Pro will select the USB connection as the main communication channel. If you want to change to the serial connection you can specify the change by pressing the Touch Cursor above *ser* in the main Local mode menu.

# Index

- BACK EAR, 11
- back side layout, 9
- BACK TAB, 11
- braille table, 21, 23
- changing of settings, 17
- communication status, 18
- config submenu, 21
- connections, 9
- current configuration, 19
- current setting, 17
- CURSOR, 10
- dot pressure, 20, 22
- Double Touch Cursor, 12, 25
- DOWN, 10
- DOWN PAD, 11
- DTC, 12
- EAST PAD, 12
- enhanced ergonomic design, 7
- Error messages, 27
- FORWARD EAR, 11
- FORWARD TAB, 11
- front panel key, 10
- front panel keys, 20, 23
- front side layout, 9
- functionality, 7
- HOME, 10
- inbound USB connector, 10
- internal braille table, 21, 23
- LEFT, 10
- LEFT PAD, 11
- load submenu, 24
- loaded set, 19
- loading settings, 17
- Local mode menu, 15, 16, 28
- main menu, 17
- main power connector, 10
- manual, 7
- MS-DOS, 7
- NORTH PAD, 11
- on/off switch, 10
- Open Standard philosophy, 7
- package, 8
- RIGHT, 10
- RIGHT PAD, 11
- satellite, 9
- satellite keypad, 10
- save submenu, 24
- saving settings, 17
- Screen reader satellite keypad, 11
- serial cable, 14
- serial connection, 14, 17, 28
- serial connector, 10
- serial number, 18
- serial port, 14
- settings menu, 19
- sleep mode, 19, 22
- SOUTH PAD, 11
- status cell, 20, 22

status submenu, 17  
switch satellites, 20, 22  
TC, 12  
test mode, 17  
test submenu, 24  
Touch Cursor, 12, 25

UP, 10  
UP PAD, 11  
USB connection, 13, 17, 28  
volume, 19, 21  
WEST PAD, 11  
Windows satellite keypad, 11